


What is the Brain Imaging Data Structure and why you should know about this!


StepuP workshop - Zürich 29/11/2024 Julius Welzel







What is the Brain Imaging Data Structure and why you should know about this!

The Past, Present, and Future of BIDS



The Past, Present, and Future of BIDS





What is BIDS?



About 20+ different EEG systems record in different file formats



Aims of BIDS




→ Makes (neuroimaging-)data FAIR

DATA STRUCTURE

- BIDS is based on simple file formats and folder structures
- Provides a "easy-to-understand" folder structure
- Provides a naming convention for files and directories
- Specifies file formats to use for modalities

METADATA

- Prevents metadata getting lost
- Some metadata is better than no metadata (80/20 rule)



GitHub repository

The Brain Imaging Data Structure

The Brain Imaging Data Structure (BIDS) is a simple and intuitive way to organize and describe data.

A Q Search

This document defines the BIDS specification, which provides many details to help implement the standard. It includes the core specification as well as many extensions to specific brain imaging modalities, and increasingly also to other kinds of data.

If BIDS is new to you, and you would like to learn more about how to adapt your own datasets to match the BIDS specification, we recommend exploring the BIDS Starter Kit. Alternatively, to get started please read the introduction to the specification.


For an overview of the BIDS ecosystem, visit the BIDS homepage. The entire specification can also be downloaded as PDF.

My favourite two BIDS principles

Requirement levels

Level of requirement for files and within files:

- 1. REQUIRED
- 2. RECOMMENDED
- 3. OPTIONAL



| participant_id | age | group |
|----------------|-----|-------|
| sub-001 | 34 | C |
| sub-002 | 12 | Р |
| sub-003 | 33 | Μ |

BIDS definitions

Dataset - A set of neuroimaging and behavioral data acquired for a purpose of a particular study. A dataset consists of data acquired from one or more subjects, possibly from multiple sessions.

Subject - A person or animal participating in the study. Used interchangeably with term **Participant**.

Session - A logical grouping of neuroimaging and behavioral data consistent across subjects.

Task - A set of structured activities performed by the participant.

Overview files:

- participant.tsv
- scans.tsv
- ...

Subject, session and task specific files

- sub-001_ses-01_task-GoNoGo.eeg
- sub-001_task-Oddball_channels.tsv

Directories and file naming structure

- Data for each subject are organized in subdirectories labeled "sub-<label>", with "<label>" representing the unique identification label of each subject
- If subjects participated in multiple sessions,
 subdirectories labeled "ses-<label>" contain session-specific data within the subject directory
- Each session subdirectory (or the subject subdirectory if no session data exists) contains subdirectories for different data types, only defined if files are present for that type
- For a data file that was collected in a given session from a given subject, the file name MUST begin with the string sub-<label>[_ses-<label>]_task-<label>

```
README.md
☐ dataset_description.json
☐ participants.json
☐ participants.tsv
☐ sub-001
☐ ses-01
☐ eeg/
☐ sub-001_ses-01_task-GoNoGo.eeg
```



```
README.md
- dataset_description.json
- participants.json
- participants.tsv
- sub-001
  — ses-01
          - eeg/
         – motion/
  — ses-02
- sub-002
```

Modality agnostic data

dataset_description.json

```
README.md
☐ dataset_description.json
☐ participants.json
☐ participants.tsv
☐ sub-001
☐ ses-01
☐ ses-02
☐ sub-002
```

```
Name": "The mother of all experiments",

"BIDSVersion": "1.4.0",

"DatasetType": "raw",

"License": "CCO",

"Authors": [ "Paul Broca", "Carl Wernicke" ],

"Funding": [ "NIH F37823MFH1" ],

"EthicsApprovals": ["HRPO (Protocol AR0928" ],

"DatasetDOI": "10.0.2.3/dfjj.10"

}
```

participants

```
README.md
☐ dataset_description.json
☐ participants.tsv
☐ participants.json
☐ sub-001
☐ ses-01
☐ ses-02
☐ sub-002
```

participants.tsv

```
participant_id age group
sub-001 34 C
sub-002 12 P
sub-003 33 M
```

participants.json

```
"age": {
    "Description": "age of the
    participant", "Units": "years"
    },

"group": {
    "Description": "assigned group",
    "Levels": {
        "C": "control",
        "M": "medication",
        "P": "placebo",
     }
}
```





Modality specific data

Files | eeg

```
sub-001
 – ses-01
    eeg/
        - sub-001 ses-01 task-<label> eeg.<extension>
        - sub-001 ses-01 task-<label> eeg.json
         sub-001 ses-01 task-<|abe|> channels.tsv
         sub-001 ses-01 task-<label> electrodes.tsv
       - sub-001_ses-01_task-<label>_electrodes.tsv
```

```
impedance
name
                          Z
CP5
        -0.77
                 -0.30
                          0.57
FC5
        -0.77
                 0.30
                          0.57
                                   12
FC1
        -0.29
                 0.31
                          0.91
                                   2
C3
        -0.59
                 0.00
                          0.81
VEOG
        n/a
                 n/a
                          n/a
                                   n/a
```

```
"TaskName": "TASKNAME",
    "SamplingFrequency": 1000,
    "SoftwareFilters": "n/a",
    "EEGChannelCount": 4,
    "EOGChannelCount": 1,
    "EEGReference": "placed on Cz",
    "PowerLineFrequency": 50
}
```

```
units
                         status
                                  status description
name
        type
CP5
        EEG
                 microV
                                  n/a
                         good
FC5
        EEG
                                  high freq noise
                 microV
                         bad
FC1
        EEG
                 microV
                         good
                                  n/a
C3
        EEG
                 microV
                                  n/a
                         good
VEOG
        EOG
                 microV
                                  n/a
                         good
```

```
{
    "EEGCoordinateSystem": "T1w",
    "EEGCoordinateUnits": "mm",
    "AnatomicalLandmarkCoordinates": {
        "LPA": [-0.067, 1.736e-09, -3.844e-09],
        "NAS": [-4.11e-09, 0.091, -4.541e-10],
        "RPA": [0.064, -6.435e-09, -4.566e-09]
    },
    "AnatomicalLandmarkCoordinateSystem": "T1w",
    "AnatomicalLandmarkCoordinateUnits": "mm",
    "IntendedFor": "sub-01_T1w.nii.gz"
}
```

Metadata | eeg.json

REQUIRED

TaskName, SamplingFrequency,
PowerlineFrequency, EEGReference,
SoftwareFilters

RECOMMENDED

TaskDescription, Manufacturer, EEGChannelCount, EOGChannelCount, ...

OPTIONAL

_

```
TaskName: Seeing stuff,
SamplingFrequency: 2400,
Manufacturer: Brain Products,
ManufacturersModelName: BrainAmp DC,
EEGChannelCount: 87,
EOGChannelCount: 2,
TriggerChannelCount: 1,
PowerLineFrequency: 50,
EEGPlacementScheme: 10 percent system,
EEGReference: single electrode placed on FCz,
SoftwareFilters: {
       Anti-aliasing filter: {
                          half-amplitude cutoff (Hz):
                           500,
                          Roll-off: 6dB/Octave
```

Metadata channels.tsv

REQUIRED

name, type, units

RECOMMENDED

_

OPTIONAL


description, sampling_frequency


| name | type | units | description | reference | status | status_description |
|---------|------|-------|--------------------------------|--------------|--------|----------------------|
| VEOG | VEOG | υV | left eye | VEOG-, VEOG+ | good | n/a |
| FDI | EMG | υV | left first dorsal interosseous | FDI-, FDI+ | good | n/a |
| Cz | EEG | υV | n/a | REF | bad | high frequency noise |
| UADC001 | MISC | n/a | envelope of audio signal | n/a | good | n/a |

Metadata electrodes.tsv

| name | X | У | Z | type | material | color |
|-------|-------|-------|--------|---------|----------|--------|
| VEOG+ | n/a | n/a | n/a | cup | Ag/AgCl | blue |
| VEOG- | n/a | n/a | n/a | cup | Ag/AgCl | white |
| FDI+ | n/a | n/a | n/a | cup | Ag/AgCl | red |
| FDI- | n/a | n/a | n/a | cup | Ag/AgCl | red |
| GND | -0.07 | 0.00 | -0.070 | clip-on | Ag/AgCl | pink |
| Cz | 0.00 | 0.07 | 0.06 | cup | Ag/AgCl | yellow |
| REF | -0.07 | -0.02 | -0.01 | cup | Ag/AgCl | grey |

| Column name | Requirement Level | Data type | Description |
|-----------------------|----------------------|-----------|---|
| name | REQUIRED | string | Name of the electrode contact point. |
| | | | Values in name MUST be unique. |
| | | | This column must appear first in the file. |
| х | REQUIRED | number | Recorded position along the x-axis. |
| | | | This column must appear second in the file. |
| у | REQUIRED | number | Recorded position along the y-axis. |
| | | | This column must appear third in the file. |
| z | REQUIRED | number or | Recorded position along the z-axis. |
| | | n/a | This column must appear fourth in the file. |
| type | RECOMMENDED | string | Type of the electrode (for example, cup, ring, clip-on, wire, needle). |
| | | | This column may appear anywhere in the file. |
| material | RECOMMENDED | string | Material of the electrode (for example, Tin, Ag/AgCl, Gold). |
| | | | This column may appear anywhere in the file. |
| impedance | RECOMMENDED | number | Impedance of the electrode, units MUST be in k0hm. |
| | | | This column may appear anywhere in the file. |
| Additional Columns | OPTIONAL | n/a | Additional columns are allowed if they are defined in the associated metadata file. |





Modality specific data

Files | motion

```
0,2634511
                                        0,092295
                                                   0,0086682
                                                              0,9305117
                                                                         0,690106
                             0.694520
                                        0,1918243
                                                                         0,885496
                                                   0,8437273
                                                              0,3975710
sub-001/
                             0,0766395
                                        0,2587211
                                                   0,5434792
                                                              0,2822837
                                                                         0,2789791
                             0,5779993
                                        0,0456141
                                                   0,0490745
                                                             0,9408899
                                                                         0,1533421
   — ses-01/
                             0.054556
                                        0,7915927
                                                   0,5871733
                                                              0,4669577
                                                                         0,9754468
                             0.966024
                                        0,1962834
                                                   0,7114406
                                                             0,3389448
                                                                         0,7194495
        -eeg/
                             0,984172
                                        0,5079461
                                                   0,1180168
                                                                         0,1753768
                                                              0,7966978
                             0,9883907
                                       0,1557346
                                                   0,8002013
                                                              0,6334882
                                                                         0,7526906
    •••
        -motion/
         — sub-001_ses-01_task-<label>_tracksys-<label>_motion.tsv
          - sub-001 ses-01 task-<label> tracksys-<label> motion.json
         — sub-001_ses-01_task-<label>_channels.tsv
```

```
{
"SamplingFrequency": 60,
"SamplingFrequencyEffective": 60.19,
"TaskName": "BIDS Motion fictive example",
"TrackingSystemName": "imu1",
"TaskDescription": "walking and talking",
"MotionChannelCount": 6,
"SubjectArtefactDescription": "n/a",
"TrackedPointsCount": 2,
"ACCELChannelCount": 3,
"GYROChannelCount": 3,
"Manufacturer": "BWSensing",
"ManufacturersModelName": "BW-imu600",
}
```

0.8098815

0,8952724

0,2326254

0,6683652

0,0480541

0,4384892

0,4886533

0,8529441

```
tracked point
                                                                                                 placement
name
                            component
                                          type
                                                                                   units
                                                                                   m/s^2
                                                                                                 right foot
imul rf acc x
                                          ACCEL
                           Х
                                                       rf
                                         ACCEL
                                                                                   m/s^2
                                                                                                 right _foot
imu1_rf_acc_y
                                                       rf
imu1_rf_acc_z
                                         ACCEL
                                                                                   m/s^2
                                                                                                 right _foot
                                                       rf
                           Z
imu1_rf_gyro_x
                                         GYRO
                                                                                   rad/s
                                                                                                 right foot
                           Х
                                                       rf
imu1_rf_gyro_y
                                         GYRO
                                                       rf
                                                                                   rad/s
                                                                                                 right foot
imu1_rf_gyro_z
                                          GYRO
                                                       rf
                                                                                   rad/s
                                                                                                 right _foot
                           Z
```

Metadata motion.json

REQUIRED

TaskName, SamplingFrequency

RECOMMENDED

TaskDescription, RecordingDuration,
MotionChannelCount, <type>ChannelCount,
SoftwareFilters, ...

OPTIONAL

Manufacuturer, RecordingSoftware, ...

```
"SamplingFrequency": 60,
"TaskName": "BIDS Motion fictive example",
"TrackingSystemName": "IMU Right Hand",
"TaskDescription": "walking and talking",
"InstitutionAddress": "Fictive address",
"InstitutionName": "Fictive Institution",
"MotionChannelCount": 18,
"Recording Duration": 4667.641106,
"SubjectArtefactDescription": "n/a",
"TrackedPointsCount": 2,
"ACCELChannelCount": 6,
"GYROChannelCount": 6,
"MAGNChannelCount": 6,
"Manufacturer": "BWSensing",
"ManufacturersModelName": "BW-IMU600",
```

Metadata channels.tsv

REQUIRED

name, component, type, tracked_point, units

RECOMMENDED

placement, reference_frame

OPTIONAL

. . .


status, status_description, sampling_frequency

| Restricted keyword list for channel type | | | | | | |
|---|--|--|--|--|--|--|
| Restricted keyword list for column type in alphabetic order. Note that upper-case is REQUIRED: | | | | | | |
| Keyword Description | | | | | | |
| Accelerometer channel, one channel for each spatial axis. Column component for the axis MUST be added to the *_channels.tsv file (x, y, or z). | | | | | | |
| Angular acceleration channel, one channel for each spatial axis. Column component for the axis MUST be added to the *_channels.tsv file (x, y, or z). | | | | | | |
| Gyrometer channel, one channel for each spatial axis. Column component for the axis MUST be added to the *_channels.tsv file (x, y, or z). | | | | | | |
| Joint angle channel between two fixed axis belonging to two bodyparts. Angle SHOULD be defined between proximal and distal bodypart in deg. | | | | | | |
| | | | | | | |

| name | component | type | tracked_point | units | reference_frame |
|-----------|-----------|-------|---------------|------------------|-----------------|
| t1_acc_x | X | ACCEL | LeftFoot | $m/s^{\Lambda}2$ | global |
| t1_acc_y | У | ACCEL | LeftFoot | $m/s^{\Lambda}2$ | global |
| t1_acc_z | Z | ACCEL | LeftFoot | $m/s^{\Lambda}2$ | global |
| t1_gyro_x | X | GYRO | LeftFoot | rad/s | global |
| t1_gyro_y | У | GYRO | LeftFoot | rad/s | global |
| t1_gyro_z | Z | GYRO | LeftFoot | rad/s | global |
| | | | | | |

Overview | motion

```
0,26345511
            0,092292015 0,008668652
                                         0,930514317
                                                        0,690193606
                                                                        0,809881135
0,694520294
            0,191824943 0,843726573
                                         0,397571025
                                                        0,88542996
                                                                        0,895276224
0,076633595
            0,258720111 0,547534792
                                         0,282283781
                                                        0,27890791
                                                                        0,232620594
            0,045616941 0,04903375
                                         0,940889749
                                                        0,153318421
                                                                        0,668360752
            0,791513927 0,587116733
                                         0,466957774
                                                        0,975446368
                                                                        0,048053341
            0.196283834 0.711044406
                                         0,338944328
                                                        0,719445195
                                                                        0,438488392
0,98417512
            0,507944361 0,1180168
                                         0,796692478
                                                        0,175376468
                                                                        0,488659533
0,98839607
            0,155737146 0,800206213
                                         0,633481382
                                                        0,752698206
                                                                        0,852943441
```



| | { | |
|---|---|--|
| | | "SamplingFrequency": 60, |
| | | "SamplingFrequencyEffective": 60.19, |
| | | "TaskName": "BIDS Motion fictive example", |
| | | "TrackingSystemName": "imu1", |
| | | "TaskDescription": "walking and talking", |
| | | "MotionChannelCount": 6, |
| | | "RecordingDuration": 10.05, |
| | | "SubjectArtefactDescription": "n/a", |
| | | "TrackedPointsCount" : 2, |
| | | "ACCELChannelCount": 3, |
| | | "GYROChannelCount": 3, |
| | | "Manufacturer": "BWSensing", |
| | | "ManufacturersModelName": "BW-imu600", |
| | } | |
| ı | | |

| onset | duration | trial_type |
|-------|----------|-------------|
| 2.45 | 0.1 | heel_strike |
| 3.81 | 0.1 | heel_strike |
| 4.95 | 0.1 | heel_strike |
| 6.11 | 0.1 | heel_strike |
| 7.24 | 0.1 | heel_strike |

| name | type | units | srate | tracked_point | component | placement |
|----------------|-------|-------|-------|---------------|-----------|-------------|
| imu1_rf_acc_x | ACCEL | m/s^2 | 222 | rf | х | right_foot |
| imu1_rf_acc_y | ACCEL | m/s^2 | 222 | rf | У | right _foot |
| imu1_rf_acc_z | ACCEL | m/s^2 | 222 | rf | Z | right _foot |
| imu1_rf_gyro_x | GYRO | rad/s | 222 | rf | Х | right _foot |
| imu1_rf_gyro_y | GYRO | rad/s | 222 | rf | У | right _foot |
| imu1_rf_gyro_z | GYRO | rad/s | 222 | rf | Z | right _foot |





Synchronize modalities

scans.tsv

```
sub-<label>/
  [ses-<label>/]
  sub-<label>[_ses-<label>]_scans.tsv
```

| Column name | Requirement Level | Data type | Description |
|-----------------------|----------------------|-----------|---|
| filename | REQUIRED | string | Relative paths to files. There MUST be exactly one row for each file. |
| | | | Values in filename MUST be unique. |
| | | | This column must appear first in the file. |
| acq_time | OPTIONAL | string | Acquisition time refers to when the first data point in each run was acquired. Furthermore, if this header is provided, the acquisition times of all files from the same recording MUST be identical. Datetime format and their anonymization are described in Units. |
| | | | This column may appear anywhere in the file. |
| Additional Columns | OPTIONAL | n/a | Additional columns are allowed. |

filename

eeg/sub-control01_task-nback.eeg
eeg/sub-control01_task-walking.eeg
motion/sub-control01_task-walking.tsv

acq_time

1877-06-15T13:45:30

1877-06-15T13:55:33

1877-06-15T13:56:27

*nback.eeg

*walking.eeg

*walking.tsv

What is the Brain Imaging Data Structure and why you should know about this!

For research


Validation and optimizing pipelines

Building a sustainable research environment

Article | Open access | Published: 09 February 2023 **EEG is better left alone**Arnaud Delorme ☑

Scientific Reports 13, Article number


Preprocessing Choices for P3 Analyses with Mobile EEG:A Systematic Literature Review and Interactive Exploration



For you

Load public data from the BIDS format

 EEGLab, FieldTrip and MNE provide importers



Converting your data to BIDS

- Reuse data from your Lab and others
- Foster collaborations



```
sub = \{'01', '02', '03', '04', '05', '06', '07', '08', '09', '10'\};
% for subject 3 the age is unknown, for subject 2 the sex is not specified
age = [11 96 nan 77 82 87 18 40 26 80];
sex = \{'f' [] 'f' 'f' 'm' 'm' 'm' 'm' 'm' \};
for subindx=1:numel(sub)
 cfg = [];
 cfg.datatype = 'eeg';
 % Load your data
 cfg.dataset = sub{subindx} + '.edf';
 % specify the output directory
 cfg.bidsroot = 'bids';
             = sub{subindx};
 cfq.sub
 % specify the information for the participants.tsv file
 % this is optional, you can also pass other pieces of info
 cfg.participants.age = age(subindx);
 cfq.participants.sex = sex{subindx};
 % provide the mnemonic and long description of the task
                    = 'changedetection';
 cfg.TaskName
 cfq.TaskDescription = 'Subjects were responding as fast as possible upon a
               change in a visually presented stimulus.';
 % these are EEG specific
 cfg.eeg.PowerLineFrequency = 50; % since recorded in the EU
 cfg.eeg.EEGReference
                           = 'M1'; % left mastoid
 data2bids(cfg);
```

end

BIDS Starterpack

Website

[https://bids-specification.readthedocs.io/en/stable/]

Validator

[https://bids-standard.github.io/bids-validator/]

BIDS Validator v1.14.6

Select a RIDS dataset to validate



Brain Imaging Data Structure v1.9.0







Brain Imaging Data Structure

The BIDS Specification

Introduction

Common principles

Modality agnostic files

Modality specific files

Derivatives

Longitudinal and multi-site studies

Glossary

BIDS Extension Proposals

Appendix

Changelog

The BIDS Starter Kit

Website

Tutorials

GitHub repository

The Brain Imaging Data Structure

The Brain Imaging Data Structure (BIDS) is a simple and intuitive way to organize and describe

This document defines the BIDS specification, which provides many details to help implement the standard. It includes the core specification as well as many extensions to specific brain imaging modalities, and increasingly also to other kinds of data.

If BIDS is new to you, and you would like to learn more about how to adapt your own datasets to match the BIDS specification, we recommend exploring the BIDS Starter Kit. Alternatively, to get started please read the introduction to the specification.

For an overview of the BIDS ecosystem, visit the BIDS homepage. The entire specification can also be downloaded as PDF.







Conclusions









Outlook: EEGManySteps





Thank you for listening carefully



Thanks to **Sein Jeung** for pushing this to completion over the past three years Thanks to all of the **BIDS Maintainers and Devs** who help to implement this Thanks to my **working group**, who have given me the freedom to work in this Thanks to all **participants** who ever provided data which is now in BIDS